



## **EYE/FACE WASH & BODY SPRAY**

### **LSE 1B**

G1 EYE/FACE WASH PUSH LEVER OPERATED – WALL MOUNTED



### **RECEPTOR BOWL**

High Visibility Yellow extra-large Cyclac receptor Bowl resistant to damage from alkalies, salt solutions, oils, most acids.

### **ATOMIZER HEADS**

Twin-cushion-flow Cyclac Atomizer Heads in Safety Green and Yellow with flip-flop dust covers that protect heads and automatically releases with water pressure.

### **ACTIVATION**

Yellow Cyclac Push plate/Foot pedal operated Stay-open Brass Ball valve with Stainless Steel Ball and Spindle.

### **STREAM CONTROL**

Steady water flow under varying pressure assured with Automatic pressure compensastors for 30-90 PSI.

### **SUPPLY**

½" BSP

Recommended pressure: 30 PSI/2 kg./ cm<sup>2</sup>

### **DRAIN**

1¼" BSP

### **NET SHIPPING WEIGHT**

3 kg