

# Exclusive hygiene mattress

## A very comfortable mattress

The Exclusive hygiene mattress is constructed using a special pressure-reducing multilayer system. A firm core is surrounded by a layer of extremely soft specially-treated Polyether foam. This construction provides the patient with the highest amount of comfort.

The Exclusive hygiene mattress is one of very few mattresses which has distinguished itself by being tested and approved, by meeting all the demands mandated by the Swedish SS 876 00 04 norm. This guarantees excellent comfort and durability, quite simply a very high-quality mattress that will last for many years.

This mattress's unique surface and multilayer construction reduce the shearing forces between the patient's skin and the underlying body tissue during movement.

The Exclusive hygiene mattress is highly suitable for preventing bedsores and related complications for all types of patients because of its excellent pressure-reducing qualities.

with  **Lentex**<sup>®</sup>



**EXCLUSIVE**  
with  **Lentex**<sup>®</sup>


Bromsövägen 3, 891 60 Örnsköldsvik  
Tel. +46 (0)860 22 28 50, www.jarven.se



CE • FREONFRI • UPPFYLLED BRANDKRÄV NT FIRE 037: A+B

# Exclusive hygiene mattress

Excellent quality

Mattress	Large Indenter 300 N	Large Indenter 600 N	Small Indenter 40 N	The method of measuring
Soft Euroflex	28.9	61.1	52.1	<p>To avoid the individual impact of different people when measuring pressure, Swerea IVF has developed pressure indenters, simulating the shape of selected parts of the human body. An objective evaluation method has been created by applying pressure equal to a medium-sized human body on the indenters and measuring with pressure mats. The method of measurement is described in a proposal to the Swedish standard with the designation SS 876 00 13</p> 
Waterlily	30.1	55.9	56.4	
Solid	30.6	46.8	53.7	
<b>Exclusive</b>	<b>35.9</b>	<b>53.0</b>	<b>59.2</b>	
Soft Comfort	37.5	53.6	60.5	
Comfort	41.9	54.5	64.5	

*All measurements in mm Hg. Low number equal low pressure.  
The number on "Large Indenter 600 N" is of special importance for heavier users.*

Please contact Järven Health Care or SWEREA/IVF for more information about the test method

## Features of the Exclusive hygiene mattress

The Exclusive hygiene mattress is a 3-layer mattress with a firm core surrounded by another layer of very soft foam (density 23+40+23 kg/m<sup>3</sup>). This guarantees that the mattress will provide for a high level of comfort for many years.

The Exclusive hygiene mattress is laminated with our unique surface layer, Lentex, which repels moisture and fluid. Lentex is a very flexible hygienic surface layer which contributes to the Exclusive hygiene mattress's excellent pressure-reducing qualities. The white colour aids early discovery of soiling on the surface of the mattress, an extra help for achieving better hygiene.

The Exclusive hygiene mattress provides protection against bedsores through its pressure-reducing qualities, suitable for all types of patients/users.

Article number	Standard sizes
5700	190-205x80x13 cm
5701	190-205x85x13 cm
5702	190-205x90x13 cm

Alternate sizes are available on special order. CE-certified. Freon free.

**Warranty:** 3 years for structural fatigue (Swedish standard B value test according to SS 876 00 11).  
5 years for the Lentex surface adhesion.

All Järven mattresses are marked with the manufacturing year and week as well as the order number. This can be found on the short side of the mattress.

### PRISMA - Your health care mattress guide

By using our PRISMA - mattress guide you will be able to choose the mattress that is right for you and your needs. For further information please contact us.

### PRISMA Your health care mattress guide.

#### COMFORT



#### HYGIENE



#### DURABILITY



#### PATIENT/USER WEIGHT



#### ENVIRONMENT



#### FIRE SAFETY REQ



**Järven Health Care**

Bromsvägen 3 • SE-891 60 Örnsköldsvik, SWEDEN  
Telephone: +46(0)660-22 28 50 • Fax: +46(0)660-824 10  
Internet: www.jarven.se • E-mail: healthcare@jarven.se